

QUICK REFERENCE GUIDE

round

see your skills in
sientra®

HSC & HSC+ breast implants, gel & saline sizers, fat transfer

Doctor _____ Surgery Date _____ Patient Name _____

Phone Number _____ Billing Zip Code _____ Email _____

Office Contact _____ Surgery Location _____



Low Projection

WTH cm	PROJ cm	VOL cc	HSC Smooth		HSC Textured		GEL SIZER	
			Catalog #	QTY	Catalog #	QTY	Catalog #	QTY
10.5	2.5	160cc	10610-160LP		20610-160LP		GS10610-160LP	
11.4	2.6	190cc	10610-190LP		20610-190LP		GS10610-190LP	
12.1	2.7	220cc	10610-220LP		20610-220LP		GS10610-220LP	
12.3	2.8	250cc	10610-250LP		20610-250LP		GS10610-250LP	
13.0	2.9	280cc	10610-280LP		20610-280LP		GS10610-280LP	
13.4	3.0	310cc	10610-310LP		20610-310LP		GS10610-310LP	
13.7	3.1	340cc	10610-340LP		20610-340LP		GS10610-340LP	
14.5	3.1	370cc	10610-370LP		20610-370LP		GS10610-370LP	
14.8	3.2	400cc	10610-400LP		20610-400LP		GS10610-400LP	
15.1	3.2	425cc	10610-425LP				GS10610-425LP	
15.4	3.2	450cc	10610-450LP		20610-450LP		GS10610-450LP	
15.6	3.3	475cc	10610-475LP				GS10610-475LP	
15.9	3.3	500cc	10610-500LP		20610-500LP		GS10610-500LP	
16.2	3.4	525cc	10610-525LP				GS10610-525LP	
16.5	3.5	550cc	10610-550LP		20610-550LP		GS10610-550LP	
16.7	3.6	575cc	10610-575LP				GS10610-575LP	
17.0	3.7	600cc	10610-600LP		20610-600LP		GS10610-600LP	
17.5	3.8	650cc	10610-650LP		20610-650LP		GS10610-650LP	
17.9	3.8	700cc	10610-700LP		20610-700LP		GS10610-700LP	



Low Plus Projection

WTH cm	PROJ cm	VOL cc	HSC Smooth		HSC+ Smooth		GEL SIZER	
			Catalog #	QTY	Catalog #	QTY	Catalog #	QTY
8.0	2.1	80cc	10610-080LPP		10710-080LPP		GS10610-080LPP	
9.0	2.4	110cc	10610-110LPP		10710-110LPP		GS10610-110LPP	
9.75	2.6	140cc	10610-140LPP		10710-140LPP		GS10610-140LPP	
10.0	2.7	155cc	10610-155LPP		10710-155LPP		GS10610-155LPP	
10.5	2.8	175cc	10610-175LPP		10710-175LPP		GS10610-175LPP	
10.75	2.9	185cc	10610-185LPP		10710-185LPP		GS10610-185LPP	
11.0	3.0	195cc	10610-195LPP		10710-195LPP		GS10610-195LPP	
11.25	3.0	210cc	10610-210LPP		10710-210LPP		GS10610-210LPP	
11.5	3.1	225cc	10610-225LPP		10710-225LPP		GS10610-225LPP	
11.75	3.2	235cc	10610-235LPP		10710-235LPP		GS10610-235LPP	
12.0	3.3	255cc	10610-255LPP		10710-255LPP		GS10610-255LPP	
12.5	3.4	280cc	10610-280LPP		10710-280LPP		GS10610-280LPP	
12.75	3.5	300cc	10610-300LPP		10710-300LPP		GS10610-300LPP	
13.0	3.6	320cc	10610-320LPP		10710-320LPP		GS10610-320LPP	
13.25	3.6	335cc	10610-335LPP		10710-335LPP		GS10610-335LPP	
13.5	3.7	355cc	10610-355LPP		10710-355LPP		GS10610-355LPP	
14.0	3.8	390cc	10610-390LPP		10710-390LPP		GS10610-390LPP	
14.5	3.9	440cc	10610-440LPP		10710-440LPP		GS10610-440LPP	



Moderate Projection

WTH cm	PROJ cm	VOL cc	HSC Smooth		HSC+ Smooth		GEL SIZER	
			Catalog #	QTY	Catalog #	QTY	Catalog #	QTY
10.3	3.2	190cc	10512-190MP		10712-190MP		GS10512-190MP	
10.7	3.2	210cc	10512-210MP		10712-210MP		GS10512-210MP	
11.1	3.2	230cc	10512-230MP		10712-230MP		GS10512-230MP	
11.4	3.3	250cc	10512-250MP		10712-250MP		GS10512-250MP	
11.8	3.4	270cc	10512-270MP		10712-270MP		GS10512-270MP	
11.9	3.5	285cc	10512-285MP		10712-285MP		GS10512-285MP	
12.0	3.6	300cc	10512-300MP		10712-300MP		GS10512-300MP	
12.3	3.7	320cc	10512-320MP		10712-320MP		GS10512-320MP	
12.4	3.8	335cc	10512-335MP		10712-335MP		GS10512-335MP	
12.5	3.9	350cc	10512-350MP		10712-350MP		GS10512-350MP	
12.7	4.0	365cc	10512-365MP		10712-365MP		GS10512-365MP	
12.8	4.1	380cc	10512-380MP		10712-380MP		GS10512-380MP	
12.9	4.2	395cc	10512-395MP		10712-395MP		GS10512-395MP	
13.0	4.2	410cc	10512-410MP		10712-410MP		GS10512-410MP	
13.2	4.2	425cc	10512-425MP		10712-425MP		GS10512-425MP	
13.4	4.2	440cc	10512-440MP		10712-440MP		GS10512-440MP	
13.7	4.3	470cc	10512-470MP		10712-470MP		GS10512-470MP	
14.0	4.3	510cc	10512-510MP		10712-510MP		GS10512-510MP	
14.8	4.3	550cc	10512-550MP		10712-550MP		GS10512-550MP	
15.4	4.4	600cc	10512-600MP		10712-600MP		GS10512-600MP	
15.8	4.5	650cc	10512-650MP		10712-650MP		GS10512-650MP	
16.1	4.7	700cc	10512-700MP		10712-700MP		GS10512-700MP	

To Order Call: 1.888.708.0808

The sale and distribution of this device is restricted to users and/or user facilities that provide information to patients about the risks and benefits of this device in the form and manner specified in the approved labeling provided by Sientra, Inc.

Sientra® is a trademark of Sientra, Inc. © 2023 Sientra, Inc. All rights reserved.
MDC-0270 R12

Moderate Plus Projection



WTH cm	PROJ cm	VOL cc	HSC Smooth		HSC Textured		HSC+ Smooth		HSC+ Textured		GEL SIZER		SALINE SIZER	
			Catalog #	QTY	Catalog #	QTY	Catalog #	QTY	Catalog #	QTY	Catalog #	QTY	Catalog #	QTY
9.3	3.7	175cc	10621-175MP		20621-175MP		10721-175MP		20721-175MP		GS10621-175MP		SZ10621-175MP	
9.7	3.9	195cc	10621-195MP		20621-195MP		10721-195MP		20721-195MP		GS10621-195MP		SZ10621-195MP	
10.0	3.9	215cc	10621-215MP		20621-215MP		10721-215MP		20721-215MP		GS10621-215MP		SZ10621-215MP	
10.5	3.9	235cc	10621-235MP		20621-235MP		10721-235MP		20721-235MP		GS10621-235MP		SZ10621-235MP	
11.0	3.9	255cc	10621-255MP		20621-255MP		10721-255MP		20721-255MP		GS10621-255MP		SZ10621-255MP	
11.4	4.0	285cc	10621-285MP		20621-285MP		10721-285MP		20721-285MP		GS10621-285MP		SZ10621-285MP	
11.7	4.0	305cc	10621-305MP		20621-305MP		10721-305MP		20721-305MP		GS10621-305MP		SZ10621-305MP	
11.8	4.4	325cc	10621-325MP		20621-325MP		10721-325MP		20721-325MP		GS10621-325MP		SZ10621-325MP	
12.1	4.4	355cc	10621-355MP		20621-355MP		10721-355MP		20721-355MP		GS10621-355MP		SZ10621-355MP	
12.5	4.5	385cc	10621-385MP		20621-385MP		10721-385MP		20721-385MP		GS10621-385MP		SZ10621-385MP	
12.9	4.6	415cc	10621-415MP		20621-415MP		10721-415MP		20721-415MP		GS10621-415MP		SZ10621-415MP	
13.2	4.6	435cc	10621-435MP		20621-435MP		10721-435MP		20721-435MP		GS10621-435MP		SZ10621-435MP	
13.3	4.6	455cc	10621-455MP		20621-455MP		10721-455MP		20721-455MP		GS10621-455MP		SZ10621-455MP	
13.6	4.8	485cc	10621-485MP		20621-485MP		10721-485MP		20721-485MP		GS10621-485MP		SZ10621-485MP	
13.7	4.9	505cc	10621-505MP		20621-505MP		10721-505MP		20721-505MP		GS10621-505MP		SZ10621-505MP	
13.8	5.0	525cc	10621-525MP		20621-525MP		10721-525MP		20721-525MP		GS10621-525MP		SZ10621-525MP	
13.9	5.1	540cc	10621-540MP		20621-540MP		10721-540MP		20721-540MP		GS10621-540MP		SZ10621-540MP	
14.1	5.2	575cc	10621-575MP		20621-575MP		10721-575MP		20721-575MP		GS10621-575MP		SZ10621-575MP	
14.3	5.3	590cc	10621-590MP		20621-590MP		10721-590MP		20721-590MP		GS10621-590MP		SZ10621-590MP	
14.6	5.3	625cc	10621-625MP		20621-625MP		10721-625MP		20721-625MP		GS10621-625MP		SZ10621-625MP	
14.8	5.4	650cc	10621-650MP		20621-650MP		10721-650MP		20721-650MP		GS10621-650MP		SZ10621-650MP	
15.1	5.5	695cc	10621-695MP		20621-695MP		10721-695MP		20721-695MP		GS10621-695MP		SZ10621-695MP	

High Projection



WTH cm	PROJ cm	VOL cc	HSC Smooth		HSC Textured		HSC+ Smooth		HSC+ Textured		GEL SIZER		SALINE SIZER	
			Catalog #	QTY	Catalog #	QTY	Catalog #	QTY	Catalog #	QTY	Catalog #	QTY	Catalog #	QTY
9.6	4.0	190cc	10621-190HP		20621-190HP		10721-190HP		20721-190HP		GS10621-190HP		SZ10621-190HP	
9.9	4.1	205cc	10621-205HP		20621-205HP		10721-205HP		20721-205HP		GS10621-205HP		SZ10621-205HP	
10.2	4.2	225cc	10621-225HP		20621-225HP		10721-225HP		20721-225HP		GS10621-225HP		SZ10621-225HP	
10.5	4.3	240cc	10621-240HP		20621-240HP		10721-240HP		20721-240HP		GS10621-240HP		SZ10621-240HP	
10.8	4.4	255cc	10621-255HP		20621-255HP		10721-255HP		20721-255HP		GS10621-255HP		SZ10621-255HP	
11.1	4.5	280cc	10621-280HP		20621-280HP		10721-280HP		20721-280HP		GS10621-280HP		SZ10621-280HP	
11.4	4.6	300cc	10621-300HP		20621-300HP		10721-300HP		20721-300HP		GS10621-300HP		SZ10621-300HP	
11.7	4.7	330cc	10621-330HP		20621-330HP		10721-330HP		20721-330HP		GS10621-330HP		SZ10621-330HP	
12.0	4.8	350cc	10621-350HP		20621-350HP		10721-350HP		20721-350HP		GS10621-350HP		SZ10621-350HP	
12.3	4.9	385cc	10621-385HP		20621-385HP		10721-385HP		20721-385HP		GS10621-385HP		SZ10621-385HP	
12.5	5.0	415cc	10621-415HP		20621-415HP		10721-415HP		20721-415HP		GS10621-415HP		SZ10621-415HP	
12.8	5.1	440cc	10621-440HP		20621-440HP		10721-440HP		20721-440HP		GS10621-440HP		SZ10621-440HP	
13.1	5.2	470cc	10621-470HP		20621-470HP		10721-470HP		20721-470HP		GS10621-470HP		SZ10621-470HP	
13.4	5.3	505cc	10621-505HP		20621-505HP		10721-505HP		20721-505HP		GS10621-505HP		SZ10621-505HP	
13.7	5.4	535cc	10621-535HP		20621-535HP		10721-535HP		20721-535HP		GS10621-535HP		SZ10621-535HP	
14.0	5.5	565cc	10621-565HP		20621-565HP		10721-565HP		20721-565HP		GS10621-565HP		SZ10621-565HP	
14.3	5.6	590cc	10621-590HP		20621-590HP		10721-590HP		20721-590HP		GS10621-590HP		SZ10621-590HP	
14.5	5.8	620cc	10621-620HP		20621-620HP		10721-620HP		20721-620HP		GS10621-620HP		SZ10621-620HP	
14.7	5.9	650cc	10621-650HP		20621-650HP		10721-650HP		20721-650HP		GS10621-650HP		SZ10621-650HP	
14.9	6.0	700cc	10621-700HP		20621-700HP		10721-700HP		20721-700HP		GS10621-700HP		SZ10621-700HP	


Xtra High Projection



WTH cm	PROJ cm	VOL cc	HSC Smooth		HSC+ Smooth		GEL SIZER	
			Catalog #	QTY	Catalog #	QTY	Catalog #	QTY
9.5	5.0	245cc	10621-245XP		10721-245XP		GS10621-245XP	
9.75	5.2	275cc	10621-275XP		10721-275XP		GS10621-275XP	
10.0	5.3	295cc	10621-295XP		10721-295XP		GS10621-295XP	
10.25	5.4	315cc	10621-315XP		10721-315XP		GS10621-315XP	
10.5	5.5	340cc	10621-340XP		10721-340XP		GS10621-340XP	
10.75	5.6	360cc	10621-360XP		10721-360XP		GS10621-360XP	
11.0	5.8	395cc	10621-395XP		10721-395XP		GS10621-395XP	
11.25	5.9	420cc	10621-420XP		10721-420XP		GS10621-420XP	
11.5	6.0	445cc	10621-445XP		10721-445XP		GS10621-445XP	
11.75	6.1	470cc	10621-470XP		10721-470XP		GS10621-470XP	
12.0	6.2	510cc	10621-510XP		10721-510XP		GS10621-510XP	

FAT TRANSFER

Viality™ Lipoaspirate Wash System



Description	Min Volume	Max Volume	Catalog #	QTY
The Viality System is a single-use device that attaches in-line to harvest fat tissue for washing, concentration and reinjection.	50 mL	1,050 mL	Viality-1400	

To Order Call: 1.888.708.0808