

Silimed Gluteal Implants





Silimed Gluteal Implant

Sientra offers the broadest range of body contouring products. Formed by a silicone elastomer envelope filled with silicone elastomer. The Oval Shaped implants feature an orientation mark for positioning. Used in reconstructive and cosmetic surgery, and provided in three different types. Supplied sterile.

Silicone Elastomer - Textured Surface



REF	VOLUME cc	DIMENSIONS cm		
		A	B	C
OVAL SHAPE				
Design by Dr. Rafael Vergara				
20922-150	225	15.3	10.3	2.3
20922-155	260	15.3	10.3	2.5
20922-160	300	15.6	10.9	2.7
20922-165	330	16.0	11.5	2.9
20922-168	385	16.3	11.9	3.2
20922-170	420	16.6	12.8	3.3



Template
REF - 90922-006

Silicone Elastomer - Textured Surface



REF	VOLUME cc	DIMENSIONS cm		
		A	B	C
Oval Shape				
Design by Dr. Abel de La Peña				
20942-135	145	13.6	10.0	2.0
20942-145	233	14.5	10.4	2.9
20942-160	275	16.0	11.6	3.0
20942-170	343	17.1	12.4	3.2
20942-175	385	17.3	12.8	3.5
20942-180	435	17.8	13.2	3.9
20942-185	540	18.4	13.8	4.1




Template
REF - 90942-007

Silicone Elastomer - Smooth Surface

REF	VOLUME cc	DIMENSIONS cm		
		A	B	C
ROUND SHAPE				
Design by Dr. Jorge Hidalgo				
932-115	250	11.6	4.0	2.6
932-120	290	12.1	4.3	2.9
932-125	340	12.5	4.5	2.9
932-130	440	13.3	5.5	4.0



Template
REF - 90932-004



GLUTEAL IMPLANT

- 1) IMPLANT
- 2) GLUTEAL MEDIUS
- 3) GLUTEAL MAXIMUS
- 4) SCIATIC NERVE



Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.