


Silimed Oval Carving Block



Silimed® Oval Carving Block Implant

The Silimed® Oval Carving Block is designed to be used in aesthetic, corrective or reconstructive surgery. Formed by a silicone envelope filled with silicone elastomer, with one end reinforced to allow the surgeon to suture and fix implant in place. Provided in a SMOOTH surface. Supplied sterile.

SMOOTH	REF	#	VOLUME cc	DIMENSIONS cm	
				A	B
	3232-022	0	3.0	1.6	2.2
	3232-025	1	5.0	1.8	2.5
	3232-030	2	13.0	2.5	3.4
	3232-040	3	23.0	3.0	4.2
	3232-045	4	29.0	3.2	4.7
	3232-050	5	46.0	4.0	5.0



Federal (USA) law restricts this device to sale by or on the order of a physician.