


Silimed Smooth Round Silicone Gel Implants



Silimed Smooth Round Silicone Gel Implants


The Silimed® Smooth Round Silicone Gel Implants are comprised of a silicone-elastomer shell, smooth-surfaced with high strength silicone gel. Provided in two profiles: HIGH and MODERATE. Supplied sterile.

HIGH PROFILE



REF	VOLUME cc	DIMENSIONS cm	
		A	B
10521-095	95	8.2	2.5
10521-115	115	8.8	2.6
10521-135	135	9.2	2.7
10521-155	155	9.7	2.9
10521-175	175	9.9	3.0
10521-195	195	10.1	3.1
10521-215	215	10.2	3.4
10521-235	235	10.7	3.4
10521-255	255	11.0	3.5
10521-275	275	11.1	3.7
10521-295	295	11.3	3.8
10521-315	315	11.4	4.0
10521-335	335	11.8	4.1
10521-355	355	11.9	4.2
10521-385	385	12.3	4.3
10521-415	415	12.7	4.4
10521-455	455	12.8	4.5
10521-495	495	13.3	4.7
10521-545	545	13.9	4.8
10521-605	605	14.3	4.8
10521-655	655	15.1	4.8
10521-695	695	15.4	4.9
10521-745	745	14.8	6.1
10521-795	795	15.1	6.2
10521-845	845	15.3	6.3
10521-895	895	15.8	6.4
10521-945	945	16.1	6.5
10521-995	995	16.4	6.6

MODERATE PROFILE



REF	VOLUME cc	DIMENSIONS cm	
		A	B
10512-080	80	8.1	2.1
10512-100	100	8.6	2.3
10512-120	120	9.2	2.4
10512-140	140	9.5	2.6
10512-160	160	10.0	2.7
10512-180	180	10.5	2.8
10512-200	200	10.8	2.9
10512-220	220	11.3	2.9
10512-240	240	11.6	3.0
10512-260	260	11.9	3.1
10512-280	280	12.2	3.2
10512-300	300	12.5	3.3
10512-330	330	12.7	3.4
10512-360	360	13.0	3.6
10512-390	390	13.2	3.8
10512-420	420	13.6	3.9
10512-450	450	14.0	4.0
10512-500	500	14.3	4.1
10512-550	550	14.8	4.3
10512-600	600	15.4	4.4
10512-650	650	15.8	4.5
10512-700	700	16.1	4.7
10512-750	750	16.0	5.0
10512-800	800	16.1	5.4
10512-850	850	16.2	5.7
10512-900	900	16.4	6.0
10512-950	950	16.6	6.1
10512-1000	1000	17.0	6.2



Caution: Investigational Device limited by Federal (USA) law to investigational use.